

Product Datasheet

Rat IP10 ELISA Kit (orb437394)

Description	Rat Interferon Gamma Induced Protein 10kDa (IP10) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P48973
Sensitivity	11.5pg/mL
Entrez	245920
Expiration Date	Please enquire.