

Product Datasheet

Rat IRAK4 ELISA Kit (orb437398)

Description	Rat Interleukin 1 Receptor Associated Kinase 4 (IRAK4) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D4A7K4
Sensitivity	5.9pg/mL
Entrez	300177
Expiration Date	Please enquire.