

Product Datasheet

Mouse IRF1 ELISA Kit (orb437400)

Description	Mouse Interferon Regulatory Factor 1 (IRF1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P15314
Sensitivity	0.058ng/mL
Entrez	16362
Expiration Date	Please enquire.