

Product Datasheet

Rat IRF4 ELISA Kit (orb437403)

Description	Rat Interferon Regulatory Factor 4 (IRF4) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	D3ZEX6
Sensitivity	0.133ng/mL
Entrez	291100
Expiration Date	Please enquire.