

Product Datasheet

Mouse IRF5 ELISA Kit (orb437405)

Description	Mouse Interferon Regulatory Factor 5 (IRF5) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P56477
Sensitivity	0.048ng/mL
Entrez	27056
Expiration Date	Please enquire.