

Product Datasheet

Mouse ITaC ELISA Kit (orb437409)

Description	Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9JHH5
Sensitivity	12.4pg/mL
Entrez	56066
Expiration Date	Please enquire.