

Product Datasheet

Human LFA1a ELISA Kit (orb437414)

Description	Human Lymphocyte Function Associated Antigen 1 Alpha (LFA1a) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P20701
Sensitivity	0.055ng/mL
Entrez	3683
Expiration Date	Please enquire.