

Product Datasheet

Mouse LTBP1 ELISA Kit (orb437420)

Description	Mouse Latent Transforming Growth Factor Beta Binding Protein 1 (LTBP1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8CG19
Sensitivity	0.053ng/mL
Entrez	268977
Expiration Date	Please enquire.