

Product Datasheet

Human LTBP2 ELISA Kit (orb437421)

Description	Human Latent Transforming Growth Factor Beta Binding Protein 2 (LTBP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q14767
Sensitivity	0.062ng/mL
Entrez	4053
Expiration Date	Please enquire.