

Product Datasheet

Rat LTb ELISA Kit (orb437422)

Description	Rat Lymphotoxin Beta (LTb) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6MG45
Sensitivity	6.2pg/mL
Entrez	361795
Expiration Date	Please enquire.