

Product Datasheet

Chicken NGF ELISA Kit (orb437436)

Description	Chicken Nerve Growth Factor (NGF) Sandwich ELISA Kit
Reactivity	Gallus
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05200
Sensitivity	10.6pg/mL
Entrez	396466
Expiration Date	Please enquire.