

Product Datasheet

Monkey NGF ELISA Kit (orb437443)

Description	Monkey Nerve Growth Factor (NGF) Sandwich ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	F6T3M7
Sensitivity	6.2pg/mL
Entrez	711681
Expiration Date	Please enquire.