

Product Datasheet

Rabbit PDGFA ELISA Kit (orb437448)

Description	Rabbit Platelet Derived Growth Factor Subunit A (PDGFA) Sandwich ELISA Kit
Reactivity	Rabbit
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P34007
Sensitivity	10.8pg/mL
Entrez	100352388
Expiration Date	Please enquire.