

Product Datasheet

Rat PDGFB ELISA Kit (orb437453)

Description	Rat Platelet Derived Growth Factor Subunit B (PDGFB) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q05028
Sensitivity	12.5pg/mL
Entrez	24628
Expiration Date	Please enquire.