

Product Datasheet

Mouse PDGFC ELISA Kit (orb437455)

Description	Mouse Platelet Derived Growth Factor C (PDGFC) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q8CI19
Sensitivity	13.5pg/mL
Entrez	54635
Expiration Date	Please enquire.