

Product Datasheet

Human PDGFD ELISA Kit (orb437456)

Description	Human Platelet Derived Growth Factor D (PDGFD) Sandwich ELISA Kit
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9GZP0
Sensitivity	4.7pg/mL
Entrez	80310
Expiration Date	Please enquire.