

Product Datasheet

Rat PDGFD ELISA Kit (orb437457)

Description	Rat Platelet Derived Growth Factor D (PDGFD) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9EQT1
Sensitivity	0.125ng/mL
Entrez	66018
Expiration Date	Please enquire.