

Product Datasheet

Human PDGFRa ELISA Kit (orb437458)

Description	Human Platelet Derived Growth Factor Receptor Alpha (PDGFRa) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P16234
Sensitivity	0.060ng/mL
Entrez	5156
Expiration Date	Please enquire.