

Product Datasheet

Mouse PDGFRL ELISA Kit (orb437462)

Description	Mouse Platelet Derived Growth Factor Receptor Like Protein (PDGFRL) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6PE55
Sensitivity	30.8pg/mL
Entrez	68797
Expiration Date	Please enquire.