

Product Datasheet

Mouse PLGF ELISA Kit (orb437465)

Description	Mouse Placenta Growth Factor (PLGF) Sandwich ELISA Kit
Reactivity	Mouse
Range	1.56-100pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49764
Sensitivity	0.56pg/mL
Entrez	18654
Expiration Date	Please enquire.