

Product Datasheet

Bovine PTN ELISA Kit (orb437467)

Description	Bovine Pleiotrophin (PTN) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21782
Sensitivity	11.4pg/mL
Entrez	280904
Expiration Date	Please enquire.