

Product Datasheet

Mouse PTN ELISA Kit (orb437469)

Description	Mouse Pleiotrophin (PTN) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P63089
Sensitivity	24pg/mL
Entrez	19242
Expiration Date	Please enquire.