

Product Datasheet

Rat PTN ELISA Kit (orb437470)

Description	Rat Pleiotrophin (PTN) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P63090
Sensitivity	12.5pg/mL
Entrez	24924
Expiration Date	Please enquire.