


Product Datasheet

Daxx antibody (orb43748)

Description	Mouse Monoclonal to Daxx.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant C-terminal part (aa 558-740) of human Daxx.
Target	DAXX
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9UER7
Entrez	1616
Expiration Date	12 months from date of receipt.



Immunocytochemistry staining of Daxx in ...