

Product Datasheet

Chicken TGFb1 ELISA Kit (orb437482)

Description	Chicken Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P09531
Sensitivity	6.2pg/mL
Entrez	100873157
Expiration Date	Please enquire.