

Product Datasheet

Equine TGFb1 ELISA Kit (orb437483)

Description	Equine Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Equine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	O19011
Sensitivity	6.3pg/mL
Entrez	100033900
Expiration Date	Please enquire.