

Product Datasheet

Guinea pig TGFb1 ELISA Kit (orb437485)

Description	Guinea pig Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Guinea pig
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Z1Y6
Sensitivity	7.6pg/mL
Entrez	100135628
Expiration Date	Please enquire.