

Product Datasheet

Rat TGFb1 ELISA Kit (orb437489)

Description	Rat Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P17246
Sensitivity	13.3pg/mL
Entrez	59086
Expiration Date	Please enquire.