

Product Datasheet

Rabbit TGFb1 ELISA Kit (orb437490)

Description	Rabbit Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, platelet-poor plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	G1SV95
Sensitivity	6.1pg/mL
Entrez	100008645
Expiration Date	Please enquire.