

Product Datasheet

Monkey TGFb1 ELISA Kit (orb437491)

Description	Monkey Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Monkey
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Sensitivity	13.4pg/mL
Entrez	574135
Expiration Date	Please enquire.