

Product Datasheet

Human TGFb2 ELISA Kit (orb437494)

Description	Human Transforming Growth Factor Beta 2 (TGFb2) competitive inhibition ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P61812
Sensitivity	8.64pg/mL
Entrez	7042
Expiration Date	Please enquire.