

Product Datasheet

Monkey TGFb2 ELISA Kit (orb437498)

Description	Monkey Transforming Growth Factor Beta 2 (TGFb2) Sandwich ELISA Kit
Reactivity	Monkey
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F6ZJX4
Sensitivity	14.7pg/mL
Entrez	707540
Expiration Date	Please enquire.