

Product Datasheet

Human TGFb3 ELISA Kit (orb437500)

Description	Human Transforming Growth Factor Beta 3 (TGFb3) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P10600
Sensitivity	34pg/mL
Entrez	7043
Expiration Date	Please enquire.