

Product Datasheet

Mouse TGF β 1 ELISA Kit (orb437506)

Description	Mouse Transforming Growth Factor Beta Induced Protein (TGF β 1) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P82198
Sensitivity	22.5pg/mL
Entrez	21810
Expiration Date	Please enquire.