

Product Datasheet

Rabbit TGF β 1 ELISA Kit (orb437508)

Description	Rabbit Transforming Growth Factor Beta Induced Protein (TGF β 1) Sandwich ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q95215
Sensitivity	0.036ng/mL
Entrez	100009230
Expiration Date	Please enquire.