



Product Datasheet

Cytokeratin 18 antibody (Biotin) (orb43751)

Description	Mouse monoclonal antibody against Cytokeratin 18 conjugated to Biotin
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, FC, ICC, IHC-P, IP, WB
Immunogen	Human breast carcinoma cell line PMC-42.
Target	Cytokeratin 18
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-6 µg/ml. Intracellular staining. Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P05783
Entrez	3875
Dilution Range	Flow cytometry: Recommended dilution: 2-6 µg/ml. Intracellular staining. Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(intracellular
S...



Western
blotting
analysis of
human
cytok...