

Product Datasheet

Human TIRAP ELISA Kit (orb437513)

Description	Human Toll Interleukin 1 Receptor Domain Containing Adaptor Protein (TIRAP) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P58753
Sensitivity	0.109ng/mL
Entrez	114609
Expiration Date	Please enquire.