

Product Datasheet

Guinea pig TNFa ELISA Kit (orb437518)

Description	Guinea pig Tumor Necrosis Factor Alpha (TNFa) Sandwich ELISA Kit
Reactivity	Guinea pig
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P51435
Sensitivity	5.5pg/mL
Entrez	100135630
Expiration Date	Please enquire.