

Product Datasheet

Human TNFaIP2 ELISA Kit (orb437520)

Description	Human Tumor Necrosis Factor Alpha Induced Protein 2 (TNFaIP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q03169
Sensitivity	0.109ng/mL
Entrez	7127
Expiration Date	Please enquire.