

## Product Datasheet

### Mouse TNFaIP6 ELISA Kit (orb437521)

|                          |                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Tumor Necrosis Factor Alpha<br>Induced Protein 6 (TNFaIP6) Sandwich ELISA<br>Kit      |
| <b>Reactivity</b>        | Mouse                                                                                       |
| <b>Range</b>             | 0.156-10ng/mL                                                                               |
| <b>Form/Appearance</b>   | Traditional                                                                                 |
| <b>Storage</b>           | All reagents should be kept according to the<br>instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                       |
| <b>Application notes</b> | serum, plasma, tissue homogenates and<br>other biological fluids.                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and<br>other biological fluids.                           |
| <b>Uniprot ID</b>        | <b>008859</b>                                                                               |
| <b>Sensitivity</b>       | 0.054ng/mL                                                                                  |
| <b>Entrez</b>            | <b>21930</b>                                                                                |
| <b>Expiration Date</b>   | Please enquire.                                                                             |