

Product Datasheet

Porcine TNFa ELISA Kit (orb437523)

Description	Porcine Tumor Necrosis Factor Alpha (TNFa) Sandwich ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P23563
Sensitivity	27.8pg/mL
Entrez	397086
Expiration Date	Please enquire.