

Product Datasheet

Monkey TNFa ELISA Kit (orb437527)

Description	Monkey Tumor Necrosis Factor Alpha (TNFa) Sandwich ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P79337
Sensitivity	6.9pg/mL
Entrez	102139631
Expiration Date	Please enquire.