

Product Datasheet

Rabbit TNFb ELISA Kit (orb437534)

Description	Rabbit Tumor Necrosis Factor Beta (TNFb) Sandwich ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P10154
Sensitivity	6.1pg/mL
Entrez	100328931
Expiration Date	Please enquire.