

Product Datasheet

Mouse TNFRSF17 ELISA Kit (orb437536)

Description	Mouse Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O88472
Sensitivity	0.101ng/mL
Entrez	21935
Expiration Date	Please enquire.