

## Product Datasheet

### Mouse TNFRSF4 ELISA Kit (orb437541)

|                          |                                                                                            |
|--------------------------|--------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Tumor Necrosis Factor Receptor Superfamily, Member 4 (TNFRSF4)<br>Sandwich ELISA Kit |
| <b>Reactivity</b>        | Mouse                                                                                      |
| <b>Range</b>             | 31.25-2000pg/mL                                                                            |
| <b>Form/Appearance</b>   | Traditional                                                                                |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.   |
| <b>Note</b>              | For research use only                                                                      |
| <b>Application notes</b> | serum, plasma, tissue homogenates and other biological fluids.                             |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                             |
| <b>Uniprot ID</b>        | <b>P47741</b>                                                                              |
| <b>Sensitivity</b>       | 13.7pg/mL                                                                                  |
| <b>Entrez</b>            | <b>22163</b>                                                                               |
| <b>Expiration Date</b>   | Please enquire.                                                                            |