

Product Datasheet

Mouse TNFRSF5 ELISA Kit (orb437542)

Description	Mouse Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P27512
Sensitivity	5.4pg/mL
Entrez	21939
Expiration Date	Please enquire.