

Product Datasheet

Rat TNFRSF5 ELISA Kit (orb437543)

Description	Rat Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q4QQW2
Sensitivity	5.2pg/mL
Entrez	171369
Expiration Date	Please enquire.