

Product Datasheet

Human TNFSF12 ELISA Kit (orb437544)

Description	Human Tumor Necrosis Factor Ligand Superfamily, Member 12 (TNFSF12) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O43508
Sensitivity	11.3pg/mL
Entrez	8742
Expiration Date	Please enquire.