

Product Datasheet

Human TNFSF13 ELISA Kit (orb437546)

Description	Human Tumor Necrosis Factor Ligand Superfamily, Member 13 (TNFSF13) Sandwich ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	075888
Sensitivity	16.64pg/mL
Entrez	8741
Expiration Date	Please enquire.