

Product Datasheet

Mouse TNFSF13 ELISA Kit (orb437547)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily, Member 13 (TNFSF13) Sandwich ELISA Kit
Reactivity	Mouse
Range	46.87-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9D777
Sensitivity	22.7pg/mL
Entrez	69583
Expiration Date	Please enquire.